

StecoGuide (coOrdination®)

The easy way to a laboratory made Surgical template





Diagnostic modell

Model analysis

Prosthetic planning

...what could be a solution?...




Situation model Prosthetic Planning

Wax up

Regarding bite relation

*...this is how it could look alike and
here I would like the implants to stay...*



Dublikated model

from wax up

Make template

Vacuum formed foil

Other template techniques possible



Template

Finish template



Template

Fill acrylic in hollow spaces

Radioopak acrylic for
planning templates

(e. g. with 10% BaSO₄)



Drill sleeve holes

Drill on desired position

Drill in desired axis

Template drills for each sleeve type

M.27.01.B300 – StecoGuide single sleeves

M.27.02.B400 – StecoGuide outer sleeves

M.27.03.B300 – StecoGuide inner sleeves



Insert drill sleeves

Press sleeve in drilled hole

No gluing!

Use pressing tool for easier handling

M.27.01.E235 – StecoGuide single sleeves

M.27.02.E350 – StecoGuide outer sleeves

M.27.03.E200 – StecoGuide inner sleeves $d_{2,0}$

M.27.03.E235 – StecoGuide inner sleeves $d_{2,35}$

Advantages

With laboratory processes
Simple communication tool
Planning template
Surgical template

Disadvantages

No radiodiagnostic basis for
planning (only prosthetic)

No implementation of 3D Data



info@dentalsil.nl



+31 (0)495 651 214



+31 (0)495 652 474



Find us on Facebook

www.dentalsil.com