














General Survey Camlog

| Product line | X-Line | | K-Line | | Z-Line |
|-----------------------------|---|---|---|---|--------|
| Insert ø 3.30 mm |  |  |  |  | |
| Distance height REF | 3.00 mm I.33.01.X300 | 4.50 mm I.33.01.X450 | 2.00 mm I.33.01.K200 | 3.50 mm I.33.01.K350 | |
| Insert ø 3.80 mm |  |  |  |  | |
| Distance height REF | 3.00 mm I.33.02.X300 | 4.50 mm I.33.02.X450 | 2.00 mm I.33.02.K200 | 3.50 mm I.33.02.K350 | |
| Insert ø 4.30 mm |  |  |  |  | |
| Distance height REF | 3.00 mm I.33.03.X300 | 4.50 mm I.33.03.X450 | 2.00 mm I.33.03.K200 | 3.50 mm I.33.03.K350 | |
| Insert ø 5.00 mm |  |  |  |  | |
| Distance height REF | 3.00 mm I.33.04.X300 | 4.50 mm I.33.04.X450 | 2.00 mm I.33.04.K200 | 3.50 mm I.33.04.K350 | |

| | | | |
|--|--|---|--|
| Denture magnet |  |  | |
| Distance height REF | 2.65 mm U.00.01.X265R | 5.00 mm U.00.01.K500 | |
| Positioning cuff/ Resilience ring |  |  |    |
| Resilience REF | 0.3 mm P.00.01.X1 | 0.3 mm P.00.01.X2 | 0.00 mm P.00.01.K1 |
| | | | 0.00 mm P.00.01.K2 |
| | | | 0.3 mm P.00.05.K1 |
| Impression post |  |  | |
| Distance height REF | 6.95 mm A.00.02.X695 | Please use the denture magnet for impression! | |
| Laboratory replica |  |  Stainless steel |  |
| Distance height REF | 9.00 mm M.00.01.X900 | 9.00 mm M.00.05.X900 | 7.50 mm M.00.01.K750 |
| Torque wrench insert |  | | |
| | We offer a suitable torque wrench insert for the most popular ratchets. | | |