

# Order form StecoGuide titanium sleeves



K.27.03.EN15/01.19

With this I order acc. to the Steco General Terms and Conditions (see [www.steco.de](http://www.steco.de)) the signed products. Please send the products to the address written in the box beside this.

Shipping costs:  
Please ask for shipping details

Customer-No. \_\_\_\_\_

Contact person \_\_\_\_\_

Date \_\_\_\_\_ Signature \_\_\_\_\_

Address (office or laboratory stamp):

Single sleeves	Size in mm	REF	PU	Price per PU	Quantity per PU
Inner-ø 2.00 mm	D 3.00 / d <b>2.00</b> / L5	M.27.01.D200L5	10		
Inner-ø 2.35 mm	D 3.00 / d <b>2.35</b> / L5	M.27.01.D235L5	10		
Inner-ø 2.35 mm	D 3.00 / d <b>2.35</b> / L10	M.27.01.D235	10		
Equipment:					
Template drill	D 3.00	M.27.01.B300	1		
Pressing tool	D 2.35	M.27.01.E235	1		

Inner sleeves	Size in mm	REF	PU	Price per PU	Quantity per PU
Inner-ø 1.50 mm	D 3.50 / d 1.50 / L6	M.27.03.D150L6	10		
Inner-ø 1.50 mm	D 3.50 / d 1.50 / L10	M.27.03.D150	10		
Inner-ø 1.60 mm	D 3.50 / d 1.60 / L6	M.27.03.D160L6	10		
Inner-ø 2.00 mm	D 3.50 / d 2.00 / L6	M.27.03.D200L6	10		
Inner-ø 2.00 mm	D 3.50 / d 2.00 / L10	M.27.03.D200	10		
Inner-ø 2.20 mm	D 3.50 / d 2.20 / L6	M.27.03.D220L6	10		
Inner-ø 2.35 mm	D 3.50 / d 2.35 / L6	M.27.03.D235L6	10		
Inner-ø 2.35 mm	D 3.50 / d 2.35 / L10	M.27.03.D235	10		
Inner-ø 2.50 mm	D 3.50 / d 2.50 / L6	M.27.03.D250L6	10		
Inner-ø 2.80 mm	D 3.50 / d 2.80 / L6	M.27.03.D280L6	10		
Inner-ø 1.16 mm with depth stop	D 3.50 / d 1.16 / L5	M.27.24.D116L5	10		
Inner-ø 1.30 mm with depth stop	D 3.50 / d 1.30 / L5	M.27.24.D130L5	10		
Inner-ø 2.00 mm with depth stop	D 3.50 / d 2.00 / L5	M.27.24.D200L5	10		
Inner-ø 2.20 mm with depth stop	D 3.50 / d 2.20 / L5	M.27.24.D220L5	10		
Inner-ø 2.35 mm with depth stop	D 3.50 / d 2.35 / L5	M.27.24.D235L5	10		
Inner-ø 1.00 mm for Endo	D 3.50 / d 1.00 / L5	M.27.28.D100L5	10		
Equipment:					
Template drill	D 3.50	M.27.03.B350	1		
Pressing tool	d 2.00	M.27.03.E200	1		
Pressing tool	d 2.35	M.27.03.E235	1		

Outer sleeves	Size in mm	REF	PU	Price per PU	Quantity per PU
Inner-ø 3.50 mm	D 4.00 / d 3.50 / L5	M.27.02.D350L5	10		
Inner-ø 3.50 mm	D 4.00 / d 3.50 / L6	M.27.02.D350	10		
Inner-ø 3.50 mm, open	D 5.00 / d 3.50 / L6	M.27.18.D350	5		
Equipment:					
Template drill	D 4.00	M.27.02.B400	1		
Pressing tool	d 3.50	M.27.02.E350	1		
Sleeve holder	L 6.00	M.27.15.H350	1		

Reference balls	Size in mm	REF	PU	Price per PU	Quantity per PU
Diameter 5.00 mm	D 5.00	M.27.09.D500	10		
Diameter 2.50 mm	D 2.50	M.27.09.D250	3		

D = Outer diameter / d = Inner diameter / L = Length

The prices are exclusive of shipping and tax.

# Order form StecoGuide titanium sleeves



K.27.03.EN15/01.19

With this I order acc. to the Steco General Terms and Conditions (see [www.steco.de](http://www.steco.de)) the signed products. Please send the products to the address written in the box beside this.

Shipping costs:  
Please ask for shipping details

Customer-No. \_\_\_\_\_

Contact person \_\_\_\_\_

Date \_\_\_\_\_ Signature \_\_\_\_\_

Address (office or laboratory stamp):

Guide sleeves for	Size in mm	REF	PU	Price per PU	Quantity per PU
Ankylos ND	D 5.50 / d 4.48 / L 4	M.27.15.D448	5		
Ankylos WD	D 5.60 / d 4.90 / L 4	M.27.15.D490	5		
Astra EV "ND"	D 5.60 / d 4.60 / L 4	M.27.15.D460	5		
Astra EV "WD"	D 6.00 / d 5.20 / L 4	M.27.15.D520	5		
Astra TX "small"	D 5.70 / d 4.70 / L 4	M.27.15.D470	5		
Astra TX "large"	D 6.70 / d 5.70 / L 4	M.27.15.D570	5		
BioHorizons Yellow	D 5.30 / d 4.30 / L 6	M.27.15.D432	5		
BioHorizons Green	D 6.00 / d 5.10 / L 6	M.27.15.D510	5		
BioHorizons Blue	D 7.00 / d 6.30 / L 6	M.27.15.D630	5		
Camlog ø 3.50 mm	D 4.50 / d 3.50 / L 3	M.27.15.D350	5		
Camlog ø 4.50 mm	D 5.50 / d 4.50 / L 3	M.27.15.D450L3	5		
Camlog ø 4.50 mm	D 5.50 / d 4.50 / L 4	M.27.15.D450	5		
MIS narrow	D 4.60 / d 4.00 / L 4	M.27.15.D400	5		
MIS	D 6.50 / d 5.50 / L 4	M.27.15.D550	5		
Nobel Biocare NP	D 4.70 / d 4.10 / L 3.5	M.27.15.D410	5		
Nobel Biocare RP	D 6.00 / d 5.00 / L 3.5	M.27.15.D500	5		
Nobel Biocare WP	D 7.00 / d 6.20 / L 3.5	M.27.15.D620	5		
Straumann Guided Pilot	D 3.80 / d 2.20 / L 6	M.27.26.D220L6	5		
Straumann ø 2.80 mm	D 3.80 / d 2.80 / L 6	M.27.15.D280	5		
Straumann ø 5.00 mm	D 6.00 / d 5.00 / L 5	M.27.15.D500L5	5		
Sweden & Martina ø 4.25 mm	D 5.00 / d 4.25 / L 5	M.27.15.D425	5		
Sweden & Martina ø 5.50 mm	D 6.35 / d 5.55 / L 5	M.27.15.D555	5		
XIVE ND	D 5.50 / d 4.48 / L 4	M.27.15.D448	5		
XIVE WD	D 6.00 / d 5.20 / L 4	M.27.15.D520	5		
Zimmer A	D 5.20 / d 4.20 / L 4	M.27.15.D420	5		
Zimmer B	D 6.30 / d 5.30 / L 4	M.27.15.D530	5		

Equipment:

Sleeve holder ø 2.80 mm	d 2.80	M.27.15.H280	1		
Sleeve holder ø 4.10 mm	d 4.10	M.27.15.H410	1		
Sleeve holder ø 4.48 mm	d 4.48	M.27.15.H448	1		
Sleeve holder ø 4.50 mm	d 4.50	M.27.15.H450	1		
Sleeve holder ø 4.90 mm	d 4.90	M.27.15.H490	1		
Sleeve holder ø 5.00 mm	d 5.00	M.27.15.H500	1		
Sleeve holder ø 5.20 mm	d 5.20	M.27.15.H520	1		
Sleeve holder ø 6.20 mm	d 6.20	M.27.15.H620	1		

Sleeves for anchor pins	Size in mm	REF	PU	Price per PU	Quantity per PU
Inner-ø 1.50 mm	D 3.50 / d 1.50 / L10	M.27.20.D150L10	10		

Thommen Medical double sleeve	Size in mm	REF	PU	Price per PU	Quantity per PU
Inner sleeve VECTODrill ø 2.0 mm	D 3.55 / d 2.02 / L6	M.27.25.D200L6	10		
Inner sleeve VECTODrill ø 2.8 mm	D 3.55 / d 2.88 / L6	M.27.25.D280L6	10		
Outer sleeve VECTODrill ø 3.5 mm	D 4.40 / d 3.55 / L6	M.27.25.D350L6	5		

D = Outer diameter / d = Inner diameter / L = Length

The prices are exclusive of shipping and tax.