

With this I order acc. to the Steco General Terms and Conditions (see www.steco.de) the signed products. Please send the products to the adress written in the box beside this.

Shipping costs:
Please ask for shipping details


Customer-No. _____

Contact person _____


Date _____ Signature _____

Address (office or laboratory stamp):

ATEC Drill 1.0 mm / Endoseal




REF O.27.28.B044.051



REF O.27.28.B044.052

StecoGuide Endo-Sleeve for ATEC Drill



REF M.27.28.D100L5

REF	PU	Price per PU	Quantity per PU
O.27.28.B044.051	1		
O.27.28.B044.052	1		
M.27.28.D100L5*	10		

The prices are exclusive of shipping and tax.

System expansion:

Further drilling steps can be guided by the StecoGuide outer and inner sleeve.
With the tube-in-tube concept you can guide different drill diameter in one template.
We offer various sleeves with diameters from 1.16 mm to 3.50 mm.
The titanium sleeves overview is located at www.steco.de/download.

